



CONWAY MEDICAL CENTER

PROCEDURE

When in hard copy form, refer to Policy Manager to validate this as the most current revision.

TITLE:	LAB-SPC-BBK-07.04 Blood Bank Sample Requirements-PRO		
ISSUED BY:	Sr Tech, BBK	REFERENCE #:	LAB-SPC-BBK-07.04-PRO
APPROVED BY:	Lab Director	EFFECTIVE DATE:	2002-09-06

SCOPE: All collectors and technical staff involved with Blood Banking.

PROCEDURE:

I. Patient Sample Requirement

- A. PREFERRED: pink stopper K2EDTA tube specifically for blood bank use.
- B. ALTERNATIVE when Pink not an option: Plain, no additive red stopper tube
- C. SHARING OF TUBES: BBK tubes are dedicated. NO sharing of tubes allowed.

II. Patient Identification:

- A. Labeled at Bedside
- B. Two independent patient identifiers—using patient armband and order.
- C. Scan patient wristband
- D. Collector name and time of collection, if on demo label

III. Units Being Held for Transfusion

- A. AABB three day rule—Units for transfusion will be held and used according to the AABB 3 day rule from the collection time.
- B. Post three (3) days: Repeat Type and Crossmatch .

NOTE: In emergency cases a review of the patient’s condition and blood bank history will be assessed by a pathologist to determine if the AABB 3 day rule can be extended

IV. Multiple Transfusions

If patients require additional units, the initial sample may be used for cross-matching up to the end of the AABB 3 day rule from the collection time. After the AABB 3 day rule has elapsed a new sample



CONWAY MEDICAL CENTER

PROCEDURE

When in hard copy form, refer to Policy Manager to validate this as the most current revision.

must be obtained.

V. ABO Non-Group Specific

In the event that a patient has received ABO non-group specific transfusions due to an emergency situation, a fresh sample must be drawn before returning the patient to their original blood group. Compatibility testing must be repeated on fresh patient sample to assess their status. The AABB 3 day rule may not apply to these patients and these situations need reviewed on an individual basis as necessary.

VI. Sample Retention

- A. Cord Blood – 7 days from sample collection
- B. Unit Segments – 7 days post transfusion
- C. Crossmatch Samples – 14 days from collection or (7 days post transfusion)
- D. Type and Screen – 7 days unless crossmatch order is added to the sample within the 72 hour guideline.
- E. ABO/Rh, Fetalscreen – 7 days from date of testing
- F. Samples from transfusion reactions will be held 7 days from the suspected reaction

VII. Pretransfusion Identification Monitoring/Process Improvement

Random screening of pre-transfusion specimens is conducted to monitor the labeling and appropriateness of pre-transfusion samples. Sample not meeting the labeling requirements are forwarded to the laboratory manager/assistant manager for review and notification and/or training if necessary.

**REFERENCES: AABB Standards for Blood Banks and Transfusion Services, 26th Edition, 2009, page 33.
AABB Technical Manual, 16th Edition, page 442.**

REFERENCE STANDARDS:

CAP GEN.40000; GEN.40016; GEN.40032; GEN.40050; GEN.40100; TRM.30800; TRM.30850

REVISION/REVIEW HISTORY:

Date	Affected Section(s)	Summary of Changes ('Reviewed' or details of change)
05/31/11jo		Reviewed/Revised. Saved electronically in lab G drive.
11/07/11lds	Format	Changed format to MCN policy manager.



CONWAY MEDICAL CENTER

PROCEDURE

When in hard copy form, refer to Policy Manager to validate this as the most current revision.

11/07/11lds	Format	Changed format to MCN policy manager.
03/28/13wm	ALL	No content changes. Newest format.
03/03/15wm	Section VII	Took out QA/QI coordinator also reviewed all
03/17/17wmllds	All	Reviewed no changes
09/14/17wm	Section II	Removed Meditech