

When in hard copy form, refer to Policy Manager to validate this as the most current revision

TITLE:	LAB-SPC-03.03 Microbiology Test Catalog-PRO		
ISSUED BY:	Administrative Lab Director	REFERENCE #:	LAB-SPC-03.03-PRO
APPROVED BY:	Lab Director	EFFECTIVE DATE:	2013-04-29

SCOPE: All potential staff than may order microbiology testing to be performed

PROCEDURE: All testing is to be ordered in Order Entry under MIC category

No revisions to this document...download of entire document when there are changes.

Dictionary Download 03/08/17

NAME	MNEMONIC
ABSCSS CULTURE COLLECTION INSTRUCTIONS: STERILE CUP; RED TOP CULTURETTE	ABSC
AFB CULTURE COLLECTION INSTRUCTIONS: STERILE CUP	AFB
AFB CULTURE BODY FLUID COLLECTION INSTRUCTIONS: STERILE CUP	AFBBF
AFB SMEAR COLLECTION INSTRUCTIONS: COLLECT IN STERILE CSF TUBE	AFBSM
ANAEROBIC CULTURE COLLECTION INSTRUCTIONS:ASPIRATES-BLUE ANAEROBIC CULTURETTE;FLUIDS-STERILE CUP	ANAER
BACTERIAL MENINGITIS SCREEN COLLECTION INSTRUCTIONS: SPINAL FLD. CALL SPEEDY PAK SEND TO GRANDSTRAND.	BACTAG
BED BUG COLLECTION INSTRUCTIONS: CALL MICRO *8174* FOR INSTRUCTIONS	BEDBUG
BETA STREP CULTURE GROUP B COLLECTION INSTRUCTIONS: COLLECT USING RED CULTURETTE	GBS
BLOOD CULTURE COLLECTION INSTRUCTIONS: PEDI BOTTLES-0.5ML,YELLOW;ADULT BOTTLES-5ML, GREEN AND ORANGE	BC

NAME**MNEMONIC**

CSF CULTURE CSF
COLLECTION INSTRUCTIONS: COLLECT IN STERILE CSF TUBE

DHEC INFLUENZA Special testing DHEC INFLU
COLLECTION INSTRUCTIONS: MINI-TIP SWAB IN EMPTY TUBE-CALL LAB FOR INSTRUCTIONS

DIPHTHERIA CULTURE DIPTH
COLLECTION INSTRUCTIONS: PRIOR CONSULTATION WITH MICRO DEPT IS REQUIRED

EAR CULTURE EAR
COLLECTION INSTRUCTIONS: RED CULTURETTE; STERILE CUP IF FLUID

EYE CULTURE EYE
COLLECTION INSTRUCTIONS: MINITIP OR RED CULTURETTE

FLOW CYTOMETRY - BAL BAL-FLOW
COLLECTION INSTRUCTIONS: ***FLOW CYTOMETRY, GIVE SPECIMEN TO CYTOLOGY***

FUNGAL CULTURE FC
COLLECTION INSTRUCTIONS: RED CULTURETTE; STERILE CONTAINER

FUNGAL CULTURE BLOOD FCBLD
COLLECTION INSTRUCTIONS: USE AEROBIC OR PEDIATRIC BOTTLE

FUNGAL CULTURE BODY FLUID FCBF
COLLECTION INSTRUCTIONS: STERILE CONTAINER

FUNGAL CULTURE CSF FCCSF
COLLECTION INSTRUCTIONS: COLLECT IN STERILE CSF TUBE

FUNGAL CULTURE SKIN/HAIR/NAIL FCS
COLLECTION INSTRUCTIONS: RED CULTURETTE; OR STERILE CONTAINER

FUNGAL CULTURE SPUTUM FCSP
COLLECTION INSTRUCTIONS: STERILE CONTAINER

NAME	MNEMONIC
FUNGAL CULTURE URINE COLLECTION INSTRUCTIONS: STERILE CUP	FCUR
FUNGAL SMEAR (KOH) COLLECTION INSTRUCTIONS: RED CULTURETTE; OR STERILE CUP	FS
GC CULTURE COLLECTION INSTRUCTIONS: RED CULTURETTE, TRANSPORT ASAP, PRIOR CONSULT REQUESTED	GC
GENITAL CULTURE COLLECTION INSTRUCTIONS: RED CULTURETTE; TRANSPORT TO LAB ASAP	GEN
GIARDIA AG COLLECTION INSTRUCTIONS: STOOL;CLEAN CONTAINER,PEA SIZE MINIMUM VOLUME/ASAP TO LAB	GIARD
GRAM STAIN COLLECTION INSTRUCTIONS: CULTURESWAB, OR STERILE CONTAINER	GS
H2O CULTURE COLLECTION INSTRUCTIONS: STERILE CUP	H2O
INDIA INK COLLECTION INSTRUCTIONS: STERILE CSF TUBE	II
MACROSCOPIC PARASITE EXAM COLLECTION INSTRUCTIONS: CALL MICRO FOR INSTRUCTIONS	MPE
MIA BLOOD CULTURE COLLECTION INSTRUCTIONS: MIA BLOOD CULTURE BOTTLES,CONTACT MICROBIOLOGY DEPT.	BCMIA
MISCELLANEOUS CULTURE COLLECTION INSTRUCTIONS: RED CULTURETTE OR STERILE CONTAINER	MISC
MRSA SURVEILLANCE COLLECTION INSTRUCTIONS: COLLECT USING WHITE TOP MRSA SWAB	MRSASURV
MRSA, DR ORDER COLLECTION INSTRUCTIONS: COLLECT USING WHITE TOP MRSA SWAB	MRSA, DR

